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Our experience : patient 1

Ex vivo Chemosensitivity

Cancer type : Recurred ovarian cancer

Previous treatment : Taxol +CDDP

	<u>IC90</u>	<u>IC50</u>	<u>Ex vivo sensitivity</u>
CDDP + Adria + Cyclopho	178	60	Resistant
CDDP + Adria	60	59	Weak Sensitive
Taxol + CDDP	47	17	Partial Sensitive
Taxol + CBDCA	26	18	Partial Sensitive
Taxotere + CBDCA	32	19	Partial Sensitive
Gemcitabine + vinorelbine	52	5	Sensitive

Recommendation

1. The combination of (Gemcitabine + Vinorelbine) followed by Carboplatin.

In the second phase of treatments, Carboplatin is recommended, because Cisplatin was not effective in the previous treatments.

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
Overall survival of patients with recurrent ovarian carcinoma calculated by Kaplan-Meier statistics.

Results for ATP-TCA group ▲ and control group ●

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E-mail : study@mail.gcrl.co.kr

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1. Kurbacher CM, Untch M et al., A randomized trial of chemotherapy directed a tumor chemosensitivity assay versus physician's choice in patients with recurrent platinum-resistant ovarian adenocarcinoma Lancet Internet publication 1997
2. Cree IA, Kurbacher CM ATP based tumour chemosensitivity testing: assisting new agent development. Anti-Cancer Drugs 1999;10:431-5
3. Konecny G, Crohns C et al., Correlation of drug response with the ATP tumor chemosensitivity assay in primary FIGOIII ovarian cancer, Gynecologic Oncology, 2000;77:258-63
4. Boabang P, Kurbacher CM et al. Anti-neoplastic activity of topotecan versus cisplatin, etoposide and paclitaxel in four squamous cell cancer cell lines of the female genital tract using an ATP-Tumor chemosensitivity assay., Anti-cancer Drugs;2000:11:843-848
5. Neale MH, Myatt NE et al., Comparison of ex vivo chemosensitivity of uveal and cutaneous melanoma, Melanoma research 2001;11;6:601-609
6. Nicolantonio F, Neale MH et al. Use of an ATP-based chemosensitivity assay to design new combination of high concentration doxorubicin with other drugs for recurrent ovarian cancer. Anti-Cancer Drugs, 2002;13:625-630
7. Neuber K Treosulfan in the treatment of metastatic melanoma: from chemosensitivity testing to clinical trials., Recent Results Cancer Res. 2003;161:159-79
8. Ewhitehouse PA, Knight LA On behalf of the Prots-mouth Colorectal Cancer Multidisciplinary team. Heterogeneity of chemosensitivity of colorectal adenocarcinoma determined by modified ex-vivo ATP tumor chemosensitivity assay Anti cancer., Drugs,2003; Jun;14(5):369-75
9. Di Nicolantonio F, Neale M et al. The chemosensitivity profile of retinoblastoma. Recent Results Cancer Res,2003;161:73-80

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