

()

()446-913

314

(031)260-9600 <http://www.gcrl.co.kr>

(study@mail.gcrl.co.kr)

: 09 - 035

2009 02 16

- 1.
- 2.

1.

IFN-	<ul style="list-style-type: none"> · 가: 80,500 · : 393 CZ393 	<ul style="list-style-type: none"> · 가: 73,710 · : 477 CZ393006 · 가 : 1,162.62 	<ul style="list-style-type: none"> · : (<u>2009-10</u>) · : 3 01
Amiodarone	<ul style="list-style-type: none"> · 가: 14,740 · : 450 CX371110 · 가 : 232.54 	<ul style="list-style-type: none"> · 가: 15,110 · : 450 C4504390 · 가 : 238.27 	<ul style="list-style-type: none"> · : · : 3 2
FISH(Trisomy 1q(1q25))	<ul style="list-style-type: none"> · 가: 210,000 · : 598 CZ967 	<ul style="list-style-type: none"> · 가: 235,000 · : 598 CZ967 	<ul style="list-style-type: none"> · : · : 3 2
FISH(Telvysion 1P)	<ul style="list-style-type: none"> · 가: 210,000 · : 598 CZ967 	<ul style="list-style-type: none"> · 가: 235,000 · : 598 CZ967 	<ul style="list-style-type: none"> · : · : 3 2
Phenol	<ul style="list-style-type: none"> · : GC · / : , /35 · : 1.2 ~ 35.6 mg/L 	<ul style="list-style-type: none"> · : GC/MS · / : /15 · : (mg/ g creatinine) · : 4.5 ~ 20.7 · () Phenol : 250.0 Benzen (10 ppm): 50 	<ul style="list-style-type: none"> · : · : 3 2
Tryptase	<ul style="list-style-type: none"> · 가: 198,000 · (154 CZ154) · : Heparin Plasma, EDTA Plasma, Serum 1.0 · / : , /14 · : 1.9 ~ 13.5 ug/L 	<ul style="list-style-type: none"> · 가: 71,900 · (154 CZ154) · : Serum 0.5 · / : /10 · : 0.0 ~ 11.4 ug/L 	<ul style="list-style-type: none"> · : · (Specialty. Lab) · : 3 2

2.

ACP(Acid Phosphatase)	: 0 ~ 6.5 : 0 ~ 5.5 U/L	: < 6.6 : < 6.5 U/L	· : Up-grade · : 03 02
Pyridinoline	: Below 25.0 : Below 32.0 (nmol/mol creatinine)	·Pyridinoline : 10 ~ 25, : 13 ~ 36 ·Deoxypyridinoline : 6, : 7 (umol/mol creatinine)	· : · :
Anti ss-DNA Ab IgM	Below 15 U/mL	< 17 U/mL	· : · :

3.

FTA-ABS IgM	, ,	~	· : · : 2 27
IgG Subclass , , ,	, ,	,	· : · : 2 27
Anti Cardio IgG / IgM	,	,	· : · : 2 28