

### ◆ 검사명 변경 항목 list (총 34항목)

Seq.	GC Labs 코드	변경 전	변경 후
1	K846	Specific IgE RE211(Rabbit urine proteins)	Specific IgE E211 (rabbit urine)
2	K845	Specific IgE RE206(Rabbit serum proteins)	Specific IgE E206 (rabbit serum proteins)
3	I604	Specific IgE RE215(Pigeon feathers)	Specific IgE E215 (pigeon feathers)
4	I603	Specific IgE RE217(Ferret epithelium)	Specific IgE E217 (ferret epithelium)
5	I161	Specific IgE RF221(Coffee)	Specific IgE F221 (coffee)
6	L032	Specific IgE RF224(Poppy seed)	Specific IgE F224 (poppy seed)
7	I209	Specific IgE RF226(Pumpkin seed)	Specific IgE F226 (pumpkin seed)
8	X826	Specific IgE RF234(Vanilla)	Specific IgE F234 (vanilla)
9	I602	Specific IgE RF247(Honey)	Specific IgE F247 (honey)
10	K398	Specific IgE RF261(Asparagus)	Specific IgE F261 (asparagus)
11	K904	Specific IgE RF262(Aubergine eggplant)	Specific IgE F262 (eggplant)
12	K002	Specific IgE RF263(Green pepper)	Specific IgE F263 (green pepper (unripe))
13	K003	Specific IgE RF270(Ginger)	Specific IgE F270 (ginger)
14	L401	Specific IgE RF278(Bay leaf)	Specific IgE F278 (bay leaf)
15	I208	Specific IgE RF279(Chili pepper)	Specific IgE F279 (chili pepper)
16	I166	Specific IgE RF281(Curry)	Specific IgE F281 (curry)
17	X820	Specific IgE RF288(Blueberry)	Specific IgE F288 (blueberry)
18	I176	Specific IgE RF289(Date)	Specific IgE F289 (date)
19	G324	Specific IgE RF293(Papaya)	Specific IgE F293 (papaya)
20	I175	Specific IgE RF302(Mandarin)	Specific IgE F302 (mandarin (tangerine))
21	K421	Specific IgE RF303(Halibut)	Specific IgE F303 (halibut)
22	K821	Specific IgE RF304(Spiny lobster(Langust))	Specific IgE F304 (lobster, spiny)
23	K095	Specific IgE RF306(Lime)	Specific IgE F306 (lime)
24	X871	Specific IgE RF312(Swordfish)	Specific IgE F312 (swordfish)
25	I702	Specific IgE RF314(Snail)	Specific IgE F314 (snail)
26	K787	Specific IgE RF315(Green bean)	Specific IgE F315 (green bean)
27	K976	Specific IgE RF319(Beetroot)	Specific IgE F319 (beetroot (beta vulgaris))
28	G126	Specific IgE RF322(Red currant)	Specific IgE F322 (red currant)
29	K122	Specific IgE RF332(Mint)	Specific IgE F332 (mint (mentha piperita))
30	L431	Specific IgE RM203(Trichosporon pullulans)	Specific IgE M203 (trichosporon pullulans)
31	G032	Specific IgE RT201 (Spruce)	Specific IgE T201 (spruce)
32	N270	Specific IgE RT214(Date tree)	Specific IgE T214 (datetree)
33	L011	Specific IgE RW204(Sunflower)	Specific IgE W204 (sunflower)
34	U616	Specific IgE RT205(Elder)	Specific IgE T205 (eldertree)