

()

()446-913

314

(031)260-9600 <http://www.gcri.co.kr>

(hklee@mail.gcri.co.kr)

: 08 - 65

2008 04 14

- 1.
- 2.

1.

		가
RPR()	<ul style="list-style-type: none"> : Latex Agglutination : Serum 1.0 mL / : /1 : Non-Reactive(Below 1.0 R.U) 	<ul style="list-style-type: none"> 가: 1,600 : 460-1, C4602 가 : 25.73 : 6 1
EGFR	<ul style="list-style-type: none"> : PCR-Sequencing : Tissue(Ethanol) or Paraffin block (10 μm Slide 5) / : ~ /10 	<ul style="list-style-type: none"> 가: 200,000 : 595, CZ626 : 4 14
FISH(HER2/neu)	<ul style="list-style-type: none"> : FISH / : ~ /5 : Paraffin block (2~3 μm Slide 2) Slide 3 + 4 μm : 4 μm Slide 	<ul style="list-style-type: none"> 가: 400,000 : 598, CZ967 : 5 1

2.

(:)

UniCAP F81(, cheddar)	, /10	~ /1	4 30
C. pneumoniae IgG/IgM	/1	, /1	4 14
PCR of Major/ Minor bcr abl	/7	, /3	4 15
Real-time PCR bcr/abl(BM, WB)	/3	, /3	
PML/RARA (BM, WB)	/3	, /3	
PML/RARA	/3	, /3	
Ureaplasma urealyticum PCR	/3	, /3	
Mycoplasma hominis PCR	/3	, /3	
Mycoplasma genitalium PCR	/3	, /3	
Haemophilus ducreyi PCR	/3	, /3	
Treponema pallidum PCR	/3	, /3	
Chlamydia pneumonia PCR	/3	, /3	
Mycoplasma pneumonia PCR	/3	, /3	
Enterovirus PCR	/3	, /3	

3.

Pneumonia urinary Ag	· / : /5	· / : ~ /1	· : Set-up · : 4 19
Legionella urinary Ag			
Histamine(Plasma)	· / : , /14 · : < 150 ng/dL	· / : , /21 · : Below 0.18 ng/mL	· : · : 4 16
Fibronectin	· : Nephelometry · / : /1 · : 25~40 mg/dL	· : RID · / : /8 · : 88~338 mg/L	· : · : 5 1
SRY gene	· 가: 30,870 · : 563, CY059006 · 가 : 497.06	· 가: 24,250 · : 595 C5956006 · 가 : 390.43	· : · : 5 1
17 -OH progesterone	· : (ng/mL) 3 day~2 month : < 9.50 3 month~11years : < 3.50 12 years ~ 20 years : < 1.80	· : · :	· : · : 4 28

4. Fungus culture

Fungus (Blood)	B4139 8,430	Fungus culture	9,180 405 B4051 (147.78)	· 가
Fungus (Body fluid)	B4139 8,430			: 가
Fungus (Gastrointestinal)	B4132 9,300			· (, ,)
Fungus (Genitourinary)	B4133 8,690			: (, ,)
Fungus (Others)	B4135 10,720			· (, ,)
Fungus (Respiratory)	B4131 12,360			: (, ,)

· : 5 1 · : set-up
· : 1 (: No growth after 1 weeks), 4 (: No fungus isolated)

5.

Histamine(Urine, 24hr urine)		
Insulin Free & Bound		
Retinol binding Protein(RBP)		4 30
Chlamydia P.Ab IgG/IgM()	C. Pneumonia IgG/IgM(EIA)	4 29

6. Pyruvate

Pyruvate ()
NaF(Glucose)